

HOW TO DRAW BLOOD



Apply the **tourniquet** on the **upper** part of the patient's **arm**.



Disinfect the arm with an **alcohol swab**.



Take a sterile **syringe**.



Remove the **needle cap**.



Align the **hole** in the **needle** with the **scale** for an easy **reading**.



Insert the needle into the **vein**.



Pull on the **plunger** to start **drawing blood**.



Release the tourniquet.



Apply a **gauze pad** onto the injection site and remove the needle.



Empty the blood in a vacutainer and **discard** the needle into a **suitable container**.

(to) apply	poser
tourniquet	garrot
upper	supérieur
arm	bras
alcohol swab	tampon imbibé d'alcool
syringe	seringue
(to) remove	retirer
(needle) cap	capuchon (de l'aiguille)
hole	trou
needle	aiguille
scale	échelle
reading	lecture
vein	veine
(to) pull on	tirer sur

plunger	piston
(to) draw blood	prélever du sang
(to) release (the tourniquet)	desserrer le garrot
gauze pad	compresse de gaze
(to) empty	vider
(to) discard	jeter
suitable	approprié
container	conteneur